


Remasoft 3.0


Paleomagnetic data browser and analyzer



Remasoft 3.0
Paleomagnetic data browser and analyzer
for Windows

Written by:
Martin Chadima
chadima@sci.muni.cz
Frantisek Hrouda
fhrouda@agico.cz

(release Mar 16, 2007)



Agico, Inc.
Jecna 29a
CZ 621 00 Brno
Czech Republic
Phone: +420 541 634 323
Fax: +420 541 634 328
agico@agico.cz
www.agico.com

User manual

Release: March 31, 2007

Remasoft 3.0 - User manual

Conditions of use & resources

Remasoft 3.0 was written by Martin Chadima and František Hrouda in 2003–2007.

Remasoft 3.0 is a freeware distributed by Agico, Inc., Brno, Czech Republic. The program is provided with the purchase of Agico's JR6/JR5 Spinner Magnetometers. Current users of Agico's instruments can download Remasoft 3.0 from Agico's website:

www.agico.com

Alternatively, current version of Remasoft 3.0 is available at:

www.gli.cas.cz/chadima/Remasoft30/

Websites contain:

- Remasoft30-Install.zip – zipped archive with installation package
- Remasoft30.exe – executable version of the program (it can be used if older version of program was previously installed in the computer)
- Remasoft30-UserManual.pdf – user manual in PDF format
- Remasoft30-Examples.zip – zipped archive containing example files

Any problems & bug can be reported, and suggestion can be done by sending e-mail to: chadima@sci.muni.cz (Martin Chadima)

The use of the program should be cited as:

Chadima, M., Hrouda, F. 2006. Remasoft 3.0 – a user-friendly paleomagnetic data browser and analyzer. *Travaux Géophysiques*, **XXVII**, 20–21.

The paper **Chadima, M., Hrouda, F. Remasoft 3.0: A user-friendly paleomagnetic data browser and analyzer** was submitted to *Computers & Geoscience* in 2007.

Acknowledgement

Remasoft 3.0 software was developed and tested in strong collaboration with Paleomagnetic laboratory, Institute of Geology, Prague, Czech Republic. All lab staff is acknowledged for their help and suggestions. Special thanks are due to the scientist-in-chief Petr Pruner for his approval of this work, patience and suggestions; and to Petr Schnabl for running a number of half-functioning versions, finding numerous errors and giving many important suggestions.

User manual conventions

- File names and directories are written in Courier, e.g., C:\Program Files\Remasoft\, Remasoft30.exe
- Program menu commands are written in *Italic*, e.g., *File* → *Open*
- Program buttons are written in SMALL CAPS, e.g., SAVE, MAKE SPECIMEN FILE(S)
- Names referring to exact expressions as appeared in user interface are underlined, e.g. Coordinates selector, Preview window
- Descriptions of program actions are indicated in bolt with diamond bullet, e.g.
◆ **Selecting coordinate system**

Table of contents

1. Introduction.....	6
2. Requirements.....	7
3. Installation.....	8
4. Program structure.....	9
4.1. Preview of acquired data (Step 1).....	10
4.2. Individual specimen processing (Step 2).....	10
4.3. Group of specimens processing (Step 3).....	11
5. User interface.....	12
5.1. Preview window.....	12
5.2. Browser window.....	16
5.2.1. Tab-strips.....	18
5.2.1.1. Overview tab-strip.....	18
5.2.1.2. Module plots tab-strip.....	20
5.2.1.3. XYZ plots tab-strip.....	22
5.2.1.4. Principal component analysis (PCA) tab-strip.....	22
5.2.1.5. Great circle tab-strip.....	24
5.2.1.6. Edit tab-strip.....	25
5.2.2. Auxiliary & Tool windows.....	27
5.2.2.1. Settings window.....	27
5.2.2.2. Geological file window.....	28
5.2.2.3. Coordinate system rotation window.....	30
5.2.2.4. Plane from 2 lineations window.....	30
5.3. Group statistic window.....	32
5.3.1. Tool windows.....	37
5.3.1.1. Angle between vectors.....	37
5.3.1.2. Magnetic inclination.....	38
References.....	39

Quick actions

Preview window

◆ Preview of acquired remanence data	14
◆ Screening through virtual specimens	14
◆ Selecting coordinate system.....	14
◆ Selecting record(s).....	14
◆ Importing geological data.....	14
◆ Importing magnetic susceptibility data.....	14
◆ Making specimen files	15
◆ Closing Preview window	15
◆ Termination of program	15

Browser window

◆ Launching Browser window	17
◆ Closing Browser window.....	17
◆ Opening specimen file	17
◆ Converting external specimen file(s).....	17
◆ Selecting coordinate system.....	18
◆ Selecting record(s).....	18
◆ Exporting the currently displayed graphics into graphic file	19
◆ Printing the currently displayed graphics on system default printer	19
◆ Multiple printing.....	19
◆ Zooming the Zijderveld diagram	24
◆ Performing great circle analysis	24
◆ Importing geological data.....	25
◆ Importing magnetic susceptibility data.....	26
◆ Deleting record	27
◆ Launching Settings window	27
◆ Saving setting	28
◆ Closing Setting window without saving.....	28
◆ Launching Geological file window	28
◆ Saving geological file	29
◆ Deleting record	29
◆ Rotating vector in various coordinate systems.....	30
◆ Calculating a plane from two lineations	30

Group statistics window

◆ Launching the Group statistics window	33
◆ Closing the Group statistics window	33
◆ Opening specimen files for the groups statistics	33
◆ Closing specimen files	34
◆ Selecting demagnetization step/components/great circles	34
◆ Subsequent selecting within displayed data	36
◆ Disregarding selected record	36
◆ Transposition of the reversely-polarized data	36
◆ Saving individual data	36
◆ Exporting the graphics into graphic file	37
◆ Printing the currently displayed graphics	37
◆ Calculating angle between two vectors	37

1. INTRODUCTION

Remasoft, version 3.0, is a user-friendly browser and basic analyzer of paleomagnetic data typically acquired during alternating-field/thermal demagnetization or magnetization procedures. Even though Remasoft 3.0 was originally written for the purpose of AGICO JR6/JR5 Spinner Magnetometers (*.jr6 / *.jra data formats, respectively) it features routines for importing various common paleomagnetic data formats (e.g. 2G binary; Randy Enkin's GSC-Pacific *.pmd data formats), and column-based data stored in the ASCII text files.

After opening/importing, the demagnetization data for each specimen are stored in separate specimen files allowing for easy viewing, editing, and analyzing. Analyses on specimen level include line-find method of principal component analysis (Kirschvink 1980) and analysis of remagnetization circles (Halls 1976). After analyses on specimen level, a number of specimen files can be combined for the purpose of group statistics. Group statistics includes calculation of mean direction (and Fisher statistics, Fisher 1953), virtual geomagnetic pole (VGP), and orientation matrix. Using orientation matrix, normal- and reversed-polarity directions can be distinguished and treated separately. Fisher statistics/VGP results can be saved in ASCII text files for the purpose of further processing.

The source code of Remasoft 3.0 is written in Microsoft Visual Basic 6.0 (with SP5) with extensive use of the Win32 Application Programmer's Interface (API) functions. Using API graphic functions, all graphics can be easily displayed on screen, directly printed or exported into portable graphic formats (Windows metafile, *.wmf, Windows bitmap, *.bmp). The vector graphics stored as Windows metafile can be easily opened and edited in a number of vector graphics processing programs, e.g Corel Draw, Adobe Illustrator.

2. REQUIREMENTS

An IBM-PC computer with Microsoft Windows operating system is required (MS Windows 95, 98, NT, 2000, XP).

Screen resolution should be 1024 × 768 or larger.

For correct reading of input files the program requires a **point to be set as decimal delimiter** (e.g. 4.52). For that purpose the program automatically sets a point as default decimal delimiter. This applies for some regional settings of the operating system when default decimal delimiter is not a point (e.g. 4,52). Unless an unexpected breakdown of the program occurs the custom decimal delimiter is automatically reset when the program is terminated.

After termination the user settings are stored in setting file. For that reason it is absolutely necessary that the current user is authorized to write into the program directory (by default: C:\Program Files\Remasoft\, unless the program was installed into a custom directory).

3. INSTALLATION


To install Remasoft 3.0 run `setup.exe` from the installation package and follow the given instructions. As soon as the installation is finished you can start the program by pressing **Start** → **Programs** → **Remasoft** → **Remasoft**.

To uninstall the program go to **Settings** → **Control Panel** → **Add or Remove Programs**, select Remasoft and click on CHANGE/REMOVE button.

♥**Tip:** Once the software was installed to your computer any future update(s) can be done by overwriting the executable version in program directory (by default: `C:\Program Files\Remasoft\Remasoft30.exe`). The current executable version can be downloaded from the above-mentioned websites.

4. PROGRAM STRUCTURE

The program works in following steps:



4.1. Preview of acquired data (Step 1)

The primary input data formats for Remasoft 3.0. are AGICO *.jr6 / *.jra magnetic remanence data files acquired during the on-line measurements on a set of specimens progressively demagnetized/magnetized in a number of steps. Two options of data storage exist: separate *.jr6 / *.jra files for each demagnetization/magnetization step or one large *.jr6 / *.jra file containing data for all specimens and all demagnetization steps.

Prior opening of *.jr6 / *.jra remanence data files several user settings should be done. The settings include the user conventions for denoting demagnetization step, user units, and user Orientation parameters.

After opening of *.jr6 / *.jra file(s), virtual specimen files are created in computer memory. The course of demagnetization/magnetization for each virtual specimen file can be monitored on computer screen. **Step 1** can be repeated after each demagnetization/magnetization step or it can be performed just once after demagnetization/magnetization experiment is terminated.

In case that *.jr6 / *.jra file(s) do(es) not contain geological data (usually inserted during on-line measurement) geological data can be imported from geological data file(s) (*.ged). The geological data file can be created using a spreadsheet-like interface built in the program or MS Excel template optionally provided with the Remasoft 3.0 distribution package.

Note: The only difference between *.jr6 and *.jra files is that the later do not contain information about Orientation parameters; during opening of *.jra files default orientation parameters are used (as stored in the user setting). Another option is the import of orientation parameters together with geological data as stored in the geological data file (*.ged).

Optionally, magnetic susceptibility data acquired during on-line measurement after each demagnetization step can be imported into the virtual specimen files. Several formats of the susceptibility files are supported (e.g. Agico formats *.su3, *.su1, *.klf, Tomasz Werner data format *.kmi).

4.2. Individual specimen processing (Step 2)

When the demagnetization/magnetization experiment is terminated **Step 1** must be performed. Virtual data files are then written to disk as individual specimen files (*.rs3). Optionally, the specimen files can be directly imported from various sources (e.g. 2G binary, Randy Enkin's GSC-Pacific *.pmd data formats, or column-based data stored in ASCII text files).

Individual specimen files can be viewed, edited, and analyzed applying the basic paleomagnetic analyses. Analyses on specimen level include the line-find method of principal component analysis (Kirschvink 1980) and the analysis of remagnetization circles (Halls 1978). The results of analyses are appended as new records into the specimen files.

Additionally, geological data and magnetic susceptibility data can be imported either manually or using the geological data files and magnetic susceptibility files, respectively. Individual records can be deleted from the specimen files.

Any graphics displayed on the screen can be directly printed or exported into the graphic output files (Windows metafile, *.wmf; or Windows bitmap, *.bmp). The vector graphics stored as Windows metafile can be easily opened and edited in a number of vector graphics processing programs, e.g. Corel Draw.

4.3. Group of specimens processing (Step 3)

After analyses on individual specimen level, group statistics can be carried out. For that purpose a number of specimen files are pooled into one large virtual data file. Records corresponding to one or more demagnetization steps or principal components are then selected from the virtual data file and used as data input for the group statistics calculations. Selected data can be saved into the Fisher statistics output file (*.rsf), original Agico remanence data file (*.jrg) or into the ASCII text file which can be used for further processing by means of other software (e.g. Magnetostratigraphy MPS.exe, a freeware written by Otakar Man, Geological Institute, Prague, Czech Republic).

Group statistics includes calculation of mean direction (Fisher statistics, Fisher 1953), virtual geomagnetic pole (VGP), and orientation matrix. Using orientation matrix, reversed-polarity directions can be distinguished and transformed to normal polarity. Results of the Fisher statistics and the orientation matrix calculations can be appended into the Fisher statistics output files (*.rsf). Later, the Fisher statistics output files can be re-opened in order to display a set of mean vectors and to perform a group statistics on the mean vector (i.e. site) level. Results of the VGP calculations can be appended into the ASCII text file which can be used for further data display or processing (e.g. VGP projections facilitated by Ppop.exe, a freeware written by Otakar Man, Geological Institute, Prague, Czech Republic).

Any graphics displayed on screen can be directly printed or exported into the graphic output files (windows metafile, *.wmf).

5. USER INTERFACE

The program consists of two main windows:

- Browser window (simplified as Preview window) – browser-like interface allowing for a rapid screening through a large number of analyzed specimens.
- Group statistics window – facilitates analyses on a group of specimens.


Both Browser (Preview) and Group statistics windows can be opened simultaneously allowing for switching between displays of individual specimen and group of specimens.

In addition, several auxiliary and tools windows are built into the program:

- Settings window
- Geological file window
- Coordinate system rotation window
- Plane from 2 lineations window
- Angle between vectors window
- Magnetic inclination window

5.1. Preview window

Preview window is used in the Preview of acquired data (Program **Step 1**).



Preview window consists of:

- (1) Headline – displays full path of currently opened file(s), [in brackets].

Remasoft 3.0 - User manual

- (2) Information bar – displays basic common information of the current specimen. These include Name of specimen, Type of demagnetization, Number of Records (number of Demagnetization steps + Principal components/Great circles, in parentheses), Specimen orientation, Dip direction (Strike) / Dip of the Bedding, Trend / Plunge of the Fold axis, and four Orientation parameters.
- (3) Virtual Specimens selector
- (4) Spherical projection plot – spherical projections of vector directions after each demagnetization/magnetization step either in the equal-angle projection (the Wulf or stereographic projection, denoted as Wulf) or in the equal-area projection (the Lambert or Schmidt projection, denoted as Lambert). Default projection can be selected in the Settings window (*Tools* → *Settings*). Solid symbols represent directions in the lower hemisphere, open symbols represent directions in the upper hemisphere. Individual points are connected with straight lines, direction of the first record (usually NRM) is represented by a crossed symbol.
- (5) Zijdeveld diagram – the two-plane Zijdeveld diagram (Zijdeveld 1967). Projection into the horizontal plane is displayed in blue, projection into the vertical plane in green. The planes of projection can be chosen in the Settings window (*Tools* → *Settings*).
- (6) Module plot – diagram of intensity of magnetization measured after each demagnetization step normalized by the maximum value of magnetization (M / M_{max}) plotted against the intensity of demagnetization. The individual points are shown as black circles connected with straight black lines. Both (M / M_{max}) and demagnetization axes are labeled automatically. Optionally, the normalized intensities of the Cartesian components (M_x, M_y, M_z) in the current coordinate system are shown in magenta, green, and blue, respectively.
- (7) Magnetic susceptibility plot/Table of demagnetization data – diagram of magnetic susceptibility as measured after each demagnetization step. The individual points are shown as black circles connected with straight black lines. Optionally, Table of demagnetization data can be displayed by double-clicking on Magnetic susceptibility plot. Another double-click on the Table restores the display of Magnetic susceptibility plot.
- (8) Coordinates selector – coordinate systems available for the current specimen are highlighted. Coordinate systems are:
 - Specimen** coordinate system (abbreviated as Spec): right-hand orthogonal coordinate system (X, Y, Z axes) where X-axis correspond to the fiducial mark drawn on specimen.
 - Geographic** or **In-situ** coordinate system (abbreviated as Geo): X-axis points to the North, Y-axis points to the East and Z-axis points downwards.
 - Tilt correction** coordinate system (abbreviated as Tilt): the same as above after rotation of the bedding to the horizontal according the strike of the bedding.
 - Full tectonic correction** coordinate system (abbreviated as Full): the same as above after rotation of the fold axis to the horizontal.
- (9) Records selector – contains a list of demagnetization steps/components stored on disk in the current specimen file. It displays a numerical order of records and demagnetization states/principal component names.
- (10) IMPORT GEOLOGICAL DATA button
- (11) IMPORT SUSCEPTIBILITY button
- (12) MAKE SPECIMEN FILE(S) button

◆ **Preview of acquired remanence data**

1. Press *File* → *Open* (or *CTRL + O*) in the main menu.
2. Select type of file in the List file of types: JR6 file(s) (* .jr6), JRA file(s) (* .jra), Geological data file(s) (* .ged), original Prague JRS file(s) (* .jrs), GSC-Pacific PMD file(s) (* .pmd), 2G binary file(s) (* .dat), 2G Hawaii file(s) (* .dat), 2G Zurich file(s) (* .txt).

Note: Opening of geological data files will be described later.

3. Select file(s) by either dragging the mouse over File names (left mouse button pressed [+CTRL for individual multiple selection, +SHIFT for continuous multiple selection]) and click on the *OPEN* button.
4. List of available virtual specimens is displayed in the Specimens selector and the plots of the first specimen are shown.

◆ **Screening through virtual specimens**

Click on the desired file in the Specimens files selector. To screen through the specimens press *Up* and *Down* keyboard keys.

◆ **Selecting coordinate system**

Click on the Coordinates selector to display data in desired coordinate system.

◆ **Selecting record(s)**

1. Click on the desired record(s). Multiple selections are facilitated by clicking and dragging the mouse or click the mouse with the left mouse button pressed [+CTRL for individual multiple selection, +SHIFT for continuous multiple selection].
2. The selected record(s) is (are) highlighted in the Records selector and is (are) displayed in red color in the plots.
3. To deselect the records click on the headline of the Records selector.

◆ **Importing geological data**

1. Click on the *IMPORT GEOLOGICAL DATA* button
2. Select one or multiple geological file(s) (* .ged) and click on the *OPEN* button.
3. Imported geological data are immediately displayed in the information bar; orientation of remanence vectors is recalculated and diagrams are redrawn.

Note: The import routine works as follows. Specimen names stored in the geological file(s) are compared with the specimen name of the current specimen file character by character from the left. This consistency-from-the-left method is particularly useful when there are more specimens of the same orientation. This is the case, for example, of several specimens cut from one cylindrical sample or several cubes cut from one a large block sample. In the geological file one stores the orientations of these large block or cylindrical core samples as measured in the field regardless of how many specimens they will yield after laboratory cutting. The orientation of the core, e.g., CR20 is used for all specimens having a name starting with CR20: CR20-1, CR20-2, CR201, CR20/1.

◆ **Importing magnetic susceptibility data**


1. Click on the *IMPORT SUSCEPTIBILITY* button
2. Select one or multiple susceptibility file(s) (Agico's formats * .su3, * .su1, * .klf, Tomasz Werner's data format * .kmi) and click on the *OPEN* button.

- Imported susceptibility data are immediately displayed in the magnetic susceptibility plot or in the susceptibility column of the Table of demagnetization data.

Note: As opposed to the import of geological data, in importing of magnetic susceptibility data a full match in the specimen name and demagnetization step is required!


◆ Making specimen files

- Click on the MAKE SPECIMEN FILE(S) button
- Save specimen files window is displayed.



- Select Drive and Directory for new files; specimen files already stored in the selected directory are displayed in the Specimen files list. The names of specimen files (shown in the Specimen file names list) are derived from the names of virtual specimen file (Virtual specimen names list) by replacing all the filename unacceptable characters (e.g. '/', '\', '!') by a hyphen ('-'); the suffix * .rs3 is added. Click on the SAVE button.

- A confirmation window appears.



New specimen files are written into the selected directory on disk. To start immediate browsing through the new files click on YES. Click on NO to close the Preview window.

Note: When the specimen files of the same names already exist they are automatically overwritten by new specimen files!

◆ Closing Preview window


Press *File* → *Close* (or *CTRL + Z*) in the main menu.

◆ Termination of program

Press *File* → *Exit* (or *CTRL + X*) in the main menu.

5.2. Browser window

Browser window (Individual specimen processing, Program Step 2) consists of:



(1)Headline – displays [a full path of the current specimen file, in brackets], and (the numerical order of the specimen file / the total number of specimen files in current directory, in parentheses).

(2)Information bar – described above in the Preview window section

(3)Drives selector

(4)Directory selector

(5)Specimen files (i.e. specimen) selector

(6)REFRESH FILE LIST button

(7)Coordinates selector – described above

(8)Records selectors – described above

(9)Tab-strip – allows for switching among six different displays of the demagnetization/magnetization data of the current specimen as opened using the Specimen files (Drives, Directory) selectors. The options are:

- Overview
- Module plots
- XYZ plots
- Line-find analysis
- Great circle
- Edit data

Remasoft 3.0 - User manual

◆ Launching Browser window

Press *File* → *Browser* (or *CTRL + B*) in the main menu.

◆ Closing Browser window

Press *File* → *Close* (or *CTRL + Z*) in the main menu.


◆ Opening specimen file

1. Launch the Browser window and select Drive and Directory. Available files are shown in the Specimen files selector.
2. Click on the desired file in the Specimen files selector.
Note: The selection of a specimen file can be done in any current Tab-strip display.
3. To screen through a large number of specimen files press *Up* and *Down* keyboard keys or *CTRL+Q*, *CTRL+A* for the previous, the next specimen file in the current directory, respectively.
4. To refresh the Specimen files selector click on the REFRESH FILE LIST button

♥**Tip:** The refresh option is useful when some specimen files are renamed, copied to, deleted from the current directory simultaneously to the operation of the program.

◆ Converting external specimen specimen file(s)

1. Press *File* → *Convert specimen file(s)* (or *CTRL + I*) in the main menu. Open multiple window appears on screen.



2. Select type of file in the List file of types: original Prague JRS file(s) (* .jrs), GSC-Pacific PMD file(s) (* .pmd), 2G binary file(s) (* .dat), 2G Hawaii file(s) (* .dat), 2G Zurich file(s) (* .txt).
3. Select file(s) by either dragging the mouse over File names (left mouse button pressed [+CTRL for individual multiple selection, +SHIFT for continuous multiple selection]) or by clicking on the SELECT ALL button. To deselect all files click on the CLEAR SELECTION button.
4. A confirmation window appears. To start immediate browsing through the new files click on YES.

Remasoft 3.0 - User manual

Note: When the specimen files of the same names already exist they files are automatically overwritten by new specimen files!

◆ Selecting coordinate system

Click on the Coordinates selector to display data in desired coordinate system.

♥ **Tip:** The selection made in Coordinates selector is effective in any Tab-strip displays of the current specimen and is also maintained when another specimen file is opened.


◆ Selecting record(s)

1. Click on the desired record(s) in the Records selector. Multiple selections are facilitated by clicking and dragging the mouse or click the mouse with the left mouse button pressed [+CTRL for individual multiple selection, +SHIFT for continuous multiple selection].
2. The selected record(s) is (are) highlighted in the Records selector and is (are) displayed in red color in all Tab-strip displays. The selection is maintained during switching among different Tab-strips options and coordinate systems.
3. To deselect the records click on the headline of the Records selector.

5.2.1. Tab-strips

5.2.1.1. Overview tab-strip

Overview tab-strip offers a quick overview of the course of demagnetization visualized in:



Remasoft 3.0 - User manual


- (1) Spherical projection plot – spherical projections of vector directions after each demagnetization step (described above in the Preview window section)
- (2) Zijderveld diagram – described above
- (3) Module plot – described above
- (4) Magnetic susceptibility plot/Table of demagnetization data – described above

◆ Exporting the currently displayed graphics into graphic file

1. Press *Graphics* → *Export* (or *CTRL + E*) in the main menu.
2. Select directory, File name, Save as type (Windows metafile, * .wmf; Windows bitmap, * .bmp) and click on **SAVE**; click on **CANCEL** to abort saving.

◆ Printing the currently displayed graphics on system default printer


1. Press *Graphics* → *Quick print* (or *CTRL + P*) in the main menu.
2. A confirmation window appears.



3. Click on **OK** to print currently displayed graphics on system default printer; click on **CANCEL** to abort printing.

◆ Multiple printing

1. Press *Graphics* → *Print* in the main menu.
2. Select a printer
Note: Selected printer will be set as Windows default printer.
3. Printing graphics using: window appears.



4. Select specimens to be printer and specify printing options.
5. Click on the PRINT button to start multiple printing. To abort printing click on the CANCEL button.


5.2.1.2. Module plots tab-strip

Module plots tab-strip shows the diagram of the intensity of magnetization measured after individual demagnetization steps, (1), normalized by the maximum value of magnetization (M / M_{max}) plotted against the intensity of demagnetization. The individual points are shown as black circles connected with straight black lines. Both (M / M_{max}) and demagnetization axes are labeled automatically. In addition to the normalized module of remanent magnetization, the normalized intensities of the Cartesian components (M_x, M_y, M_z) in currently selected coordinate system are shown as magenta squares, green triangles, and blue circles, respectively. Display of the Cartesian components of magnetization can be useful, e.g., when performing the thermal demagnetization of three orthogonal components of different coercivity (Lowrie 1990).

Optionally, the magnetic susceptibility diagram as measured after each demagnetization step can be displayed as open circles connected with gray line.

Any of above-mentioned plots can be switched on/off using the check boxes in the upper right corner of the diagram.


Absolute values corresponding to current mouse position are displayed in the black status bar in the lower part of the diagram (hit the left mouse button to enlarge the cursor cross hair, (2)).



Various demagnetization values and ratios are shown below the demagnetization diagram, (3):

- M_0 – absolute value of magnetization of the first demagnetization record (usually NRM)
- M_{max} – maximum intensity of magnetization
- M_{min} – minimum intensity of magnetization
- M_x – intensity of magnetization of the last demagnetization step (x denote the step)
- M_{max} / M_0 – maximum value of magnetization normalized by M_0 (displayed in red when $M_{max} / M_0 > 1$)
- M_{min} / M_0 – minimum intensity of magnetization normalized by M_0
- M_x / M_0 – intensity of magnetization of the last demagnetization step normalized by M_0
- $M(50\%)$ – demagnetization value when intensity of magnetization decrease to 50 % of the M_0 value (usually known as the Median destructive field in AF demagnetization)
- $M(10\%)$ – demagnetization value when intensity of magnetization decrease to 10 % of the M_0 value (most of the initial magnetization was removed; in the case of thermal demagnetization it can serve as a rough estimate of the unblocking temperature)
- $M_x(0\%)$ – demagnetization value when intensity of magnetization in the direction of X-axis is zero (useful in the back-field isothermal demagnetization parallel to the X-axis of specimen where $M_x(0\%)$ represents the coercivity of remanence, H_{cr})
- K_0 – absolute value of magnetic susceptibility of the first demagnetization record (usually magnetic susceptibility in natural state)
- K_{max} – maximum intensity of magnetic susceptibility
- K_{max} / K_0 – maximum value of magnetic susceptibility normalized by K_0


5.2.1.3. XYZ plots tab-strip



XYZ plots tab-strip shows the Cartesian components of remanent magnetization in the current coordinate system projected into three orthogonal planes. The plots are denoted according to the direction of view: FRONT VIEW (projection into Y-Z plane, (1)), SIDE VIEW (projection into X-Z plane, (2)), and TOP VIEW (projection into X-Y plane, (3)). In addition, the projection into arbitrarily oriented vertical plane is shown, (4). The orientation of projection plane can be set by moving a horizontal scroll bar below the diagram ((5), movement by 5 degrees when clicked in the scroll bar area, by 1 degree when clicked on a scroll arrow). Direction of view is displayed in degrees above the diagram as well as by the yellow arrow in the TOP VIEW diagram.

5.2.1.4. Principal component analysis (PCA) tab-strip

PCA tab-strip is designed to perform the principal component analysis in the two-plane Zijderveld diagram, (1). Projection into the horizontal plane is displayed in blue, projection into the vertical plane in green. The planes of projection can be chosen in the Settings window (*Tools* → *Settings*). After inputting the upper and lower limits of demagnetization interval the straight line is fitted using the least squares method.



◆ Performing principal component analysis

1. Input limits either manually in the Limits input boxes (2) or click on desired steps in a pair of lists of demagnetization steps.
Note: Input boxes and lists are interchangeable; after selection the program automatically sets the upper and lower limits.
2. Click on Include ORIGIN (3) to include the origin of the Zijderfeld plot to a linear segment. Click on Anchor to last point (3) to anchor a linear segment to last selected point in the Zijderfeld plot (simultaneous selection of Include ORIGIN will anchor a linear segment to the origin).
3. Click on the CALCULATE button (4) to calculate a linear segment (after calculation the button caption changes from CALCULATE to REVERSE).
4. Result of the analysis is displayed as Declination (Dec), Inclination (Inc), and the Maximum angular deviation (MAD) in the Principal component frame (5) and as orange straight line in the Zijderfeld diagram.
Note: The results appear on the screen only, nothing is stored into the specimen file yet; new selection of limits clears the Principal component frame and projection of principal component into the Zijderfeld diagram.
 Additionally, the principal component orientation and the circle corresponding to the maximum angular deviation are displayed in the spherical projection of the Overview tab-strip in orange color.
5. Click on the REVERSE button (4) to obtain the reversely -polarized component.
6. To save the principal component into the specimen file insert the Name of component and optionally a Note (6) and click on the APPEND TO FILE button (7).

7. New record is added to the Records selector. The name of record is composed of the principal component prefix (as set in the Settings window) plus the name input by the user.

Additionally, the orientation and circle corresponding to the maximum angular deviation of saved principal component(s) are displayed in the Overview tab-strip spherical projection as magenta symbol(s).

8. To delete the saved record(s) go to the Edit tab-strip.

♥**Tip:** Selections made by clicking on the lists of demagnetization steps are visualized in the Records selector and thus effective in all tab-strip options. Similarly, the selection made in the Records selector is effective in the Line-find analysis option. Selection of the limits for the purpose of the Line-find analysis can be, thus, made in any currently active tab-strip option. This may be particularly helpful when the user needs to verify the selection of limits in, e.g., spherical projection, module or susceptibility plots of demagnetization data.

◆**Zooming the Zijdeveld diagram**

1. Drag a rectangular zoom region using the mouse with the left button pressed.
2. To restore a full view click anywhere in the Zijdeveld plot.

5.2.1.5. Great circle tab-strip

Great circle tab-strip is designed for displaying vector directions in a spherical projection and for fitting of great circles within a selected interval of demagnetization steps. Vector directions measured after individual demagnetization steps are plotted (in the current coordinate system) either in the equal-angle projection (the Wulf or stereographic projection, denoted as Wulf) or in the equal-area projection (the Lambert or Schmidt projection, denote as Lambert) (1). Default projection can be selected in the Settings window (*Tools* → *Settings*). Solid symbols represent directions in the lower hemisphere, open symbols represent directions in the upper hemisphere. Individual points are connected with straight lines, direction of the first record (usually NRM) is represented by a crossed symbol.

◆**Performing great circle analysis**

1. Input limits either manually in the Limits input boxes or click on desired steps in a pair of lists of demagnetization steps (2).

Note: The input boxes and lists are interchangeable; after selection the program automatically sets the upper and lower limits.

2. By default all selected points are weighted by intensity for the purpose of the great circle calculation. Click on Normalized (3) check box to give the point on the great circle equal weight. Click on Fisher mean (3) to calculate the mean vector of selected points.


3. Click on the CALCULATE button (4) to calculate a great circle / mean vector.

3. Results of the analysis are displayed as Declination (denoted as Dec), Inclination (Inc) of the pole to the great circle, and Maximum angular deviation (MAD) in the Pole to great circle frame (5) below the diagram and as orange circle and square corresponding to the projection of the great circle and pole to it, respectively.

Note: The results appear on the screen only, nothing is stored into the specimen file yet; new selection of demagnetization interval clears the Pole to great circle frame and orange symbols in the spherical projection diagram.

4. To save the pole to the great circle into the specimen file insert the Name of the great circle and optionally a Note (6), and click on the APPEND TO FILE button (7).

5. New record is added to the Records selector. The name of record is composed of the great circle prefix (as set in the Settings window) plus the name input by the user. Additionally, saved great circle(s) and pole(s) to great circle(s) are displayed in the Overview tab-strip spherical projection in magenta symbol(s).
6. To delete the saved record(s) go to the Edit tab-strip.



♥**Tip:** Selections made by clicking on the lists of demagnetization steps are visualized in the Records selector and thus effective in all tab-strip options. Similarly, the selection made in the Records selector is effective in the Great circle option. The selection of the limits for the great circle analysis can be, thus, made in any currently active tab-strip option.

5.2.1.6. Edit tab-strip

Edit tab-strip enables editing of specimen file headline, importing geological and magnetic susceptibility data into specimen file and deleting individual records from specimen file. Geological data, (1), (i.e. Orientation of specimen, Bedding and Fold axis, Orientation parameters, Site name, Latitude, Longitude, Height, Rock, Age, Formation) can be imported either manually or from geological file(s). Magnetic susceptibility measured after each demagnetization step can be imported either manually or from susceptibility file(s).

◆ Importing geological data


1. Click on the IMPORT GEOLOGICAL DATA button, (2).

2. Select one or multiple geological file(s) (* .ged) and click on the OPEN button.

3. Imported geological data are immediately displayed in headline text boxes.

Note: The import routine works as follows. Specimen names stored in the geological file(s) are compared with the specimen name of the current specimen file (displayed in the Name text box as well as in the Information bar) character by character from the left. This consistency-from-the-left method is particularly useful when there are more specimens of the same orientation. This is the case, for example, of several specimens cut from one cylindrical sample or several cubes cut from one large block sample. In the geological file one stores the orientations of these large block or cylindrical core samples as measured in the field regardless of how many specimens they yielded after laboratory cutting. The orientation of the core, e.g., CR20 is used for all specimens having a name starting with CR20: CR20-1, CR20-2, CR201, CR20/1.

4. Click on the SAVE button, (6), for recalculation of the remanence vectors and for saving new data into the specimen file.



◆ Importing magnetic susceptibility data

1. Click on the IMPORT SUSCEPTIBILITY button, (4).

2. Select one or multiple susceptibility file(s) (Agico's formats *.su3, *.su1, *.klf, Tomasz Werner's data format *.kmi) and click on the OPEN button.

3. Imported susceptibility data are immediately displayed in the susceptibility column of the data grid.

Note: As opposed to the import of geological data, in the process of importing magnetic susceptibility data a full match in the specimen name and demagnetization step is required!

4. To save the changes in the specimen file click on the SAVE button, (6).

◆ **Deleting record**

1. Click on left column of the data grid (2) containing the records number to select the record.
2. Click on the DELETE RECORD button, (5), to delete the selected record from data grid.
Note: No changes in the specimen file are effective unless the file is saved by clicking on the Save button!
3. To save the changes in the specimen file click on the SAVE button, (6).


5.2.2. Auxiliary & Tool windows

5.2.2.1. Settings window

◆ **Launching Settings window**

Press *Auxiliary* → *Setting* (or *F12*) in the main menu.

Settings window includes:



Spherical projection frame – selection of default spherical projection to be used in the Preview, Overview windows, and Great circle tab-strip. The options are: the equal-angle projection (the Wulf or stereographic projection, denoted as Wulf) or in the equal-area projection (the Lambert or Schmidt projection, denoted as Lambert).

Zijderveld diagram frame – selection of the appearance of vector component diagram to be shown in the Preview window and Overview, Line-find analysis tab-strips. Selection can be done as follows:

1. Select a Vertical plane, either oriented north-south, denoted as N-S plane; or oriented west-east, W-E plane.
2. Select the orientation of the Horizontal plane, either north in the top position (N on top) or east on the top position (E on top).

Symbols and units frame – used in opening *.jr6 / *.jra data files or in importing external data files.

- NRM symbol – denote symbol in the State column indicating Natural remanent magnetization record(s), e.g., 'NRM', 'NS', '0', 'A0', '0mT'...etc. There are two optional boxes to set NRM symbols for distinguishing between alternating-field (AF demag) or thermal (Th.demag) demagnetizations (blank by default).

Remasoft 3.0 - User manual

- Prefix/suffix – set prefix/suffix to distinguish between alternating-field (AF demag) or thermal (Th.demag) records, e.g. 'A10' – AF demagnetization using a field of 10 [Units], T500 – thermal demagnetization using temperature 500 [Units]. Leave blank when only numbers are used, e.g. '10', '500', ...etc.
- Principal component/Great circle prefixes – set prefixes for new principal component(s) or great circle(s) names as appended to the Specimen file(s) during principal component/great circle analyses.
- Units – set units for AF or thermal demagnetization. Units are stored in Specimen files and are used for labeling horizontal axes in Module plots and in the list of Records.

Default orientation parameters frame – set default orientation parameters describing specimen orientation scheme.

Note: Default orientation parameters are used in opening *.jra file or importing external data file; *.jr6 files already contain orientation parameters.

Default site position frame – default position of sampling site to be used in the Virtual Geomagnetic Pole calculation (Group statistics window).

- Latitude – geographic latitude of the sampling site.
- Longitude – geographic longitude of the sampling site

Note: Sign conventions for geographic locations: Latitudes increase from -90° at south geographic pole to 0° at equator and to $+90^\circ$ at the north geographic pole. Longitudes east of the Greenwich meridian are positive, while westerly longitudes are negative.

Show tips – indicates whether quick tips are shown when the mouse is paused over some user interface features (i.e. diagram, selector, ...etc).

◆ Saving setting

Click on the SAVE button, setting are written into the setting file.

◆ Closing Setting window without saving

Click on the CLOSE button, window is closed without saving.

5.2.2.2. Geological file window

Geological file window is designed for creating/opening of geological data file(s).

◆ Launching Geological file window

Press *File* → *New* → *Geological file* (or *CTRL + N*) or *File* → *Open* (or *CTRL + O*) and select geological file name(s).

Remasoft 3.0 - User manual

The screenshot shows the 'Geological file(s) - [D:\antarktida\]' window. It is divided into three main sections at the top: 'Sampling site', 'Rock characteristics', and 'Orientation parameters'. Below these are input fields for site name, latitude, longitude, rock type, age, formation, and four orientation parameters (P1-P4). A large table below contains 16 rows of specimen data with columns for specimen ID, azimuth, plunge, and various orientation codes (Code1, Foli1D, Foli1I, Line1D, Line1I, Code2, Foli2D, Foli2I, Line2D, Line2I). At the bottom, there are buttons for 'Delete record', 'Save as', and 'Close'.

#	Specimen	Azimuth	Plunge	Code1	Foli1D	Foli1I	Line1D	Line1I	Code2	Foli2D	Foli2I	Line2D	Line2I
1	JR37/1A	60	10	00	0	0	0	0	00	0	0	60	60
2	JR37/2A	60	10	00	0	0	0	0	00	0	0	60	60
3	JR51/1A	60	10	00	0	0	0	0	00	0	0	60	60
4	JR51/2A	60	10	00	0	0	0	0	00	0	0	60	60
5	JR52/1A	60	10	00	0	0	0	0	00	0	0	60	60
6	JR52/2A	50	10	00	0	0	0	0	00	0	0	50	50
7	JR64/1A	50	10	00	0	0	0	0	00	0	0	50	50
8	JR64/2A	50	10	00	0	0	0	0	00	0	0	50	50
9	JR84/1A	50	10	00	0	0	0	0	00	0	0	50	50
10	JR84/2A	50	10	00	0	0	0	0	00	0	0	50	50
11	JR109/1A	50	10	00	0	0	0	0	00	0	0	50	50
12	JR109/2A	50	10	00	0	0	0	0	00	0	0	50	50
13	JR109/3A	50	10	00	0	0	0	0	00	0	0	50	50
14	JR109/4A	50	10	00	0	0	0	0	00	0	0	50	50
15	Z10/1A	50	10	00	0	0	0	0	00	0	0	50	50
16	Z10/2A	50	10	00	0	0	0	0	00	0	0	50	50
17													
18													
19													
20													
21													
22													
23													
24													

General information about sampling site are displayed in three frames:

- Sampling site information: Site name, Latitude, Longitude
- Rock characteristics: Rock name, rock Age, abbreviated Formation (Frm.)
- Orientation parameters

Orientation of individual samples are displayed in the table. Each records contains:

Specimen name, Azimuth of specimen, Plunge of specimen, and orientations of two mesoscopic foliations and two mesoscopic lineations. Each mesoscopic feature is label using a two-letter Code1 and Code2, for the 1st and 2nd foliation, respectively. Dip direction (Strike) and Dip of the foliation is denoted as Foli(1/2)D, and Foli(1/2)I, respectively. Trend and Plunge of the lineation is denoted as Line(1/2)D, and Line(1/2)I, respectively.

◆ Saving geological file

Click on the SAVE AS button, select file name, and click on the SAVE button

◆ Deleting record

1. Click on left column of the data grid containing the records number to select the record.
2. Click on the DELETE RECORD button to delete the selected record from data grid.


Note: No changes in the specimen file are effective unless the file is saved by clicking on the SAVE AS button!

5.2.2.3. Coordinate system rotation window

Coordinate system rotation window facilitates rotation of the unit vector from one coordinate system to another using orientation angle and AGICO orientation parameters. This may be particularly useful in understanding the basic idea of the orientation parameters.

◆ Rotating vector in various coordinate systems

1. Press *Tools* → *Coordinate system rotation* in the main menu.
2. Coordinate system rotation window appears.




3. Set the custom Orientation parameters.
4. Select the prior-rotation Coordinates.
5. Manually insert Vector orientation as Declination and Inclination into green text boxes.
6. Click on the DISPLAY button, vector orientation is displayed as a green point in the spherical projection diagram.
7. Manually insert Specimen orientation and/or Bedding orientation and/or Fold axis orientation.
8. Select the desired post-rotational Coordinates.
9. Rotation from the prior-rotation to post-rotation coordinate system is performed; orientation of the rotated vector is displayed in green Vector orientation text boxes and as a green point in the spherical projection diagram.
10. To clear input boxes click on the CLEAR button.

5.2.2.4. Plane from 2 lineations window

This tool window is useful, e.g., in creating a geological data file when a plane is expressed as two lineations.

◆ Calculating a plane from two lineations


1. Press *Tools* → *Plane from 2 lineations* in the main menu.
2. Plane from 2 lineations window appears.



3. Manually insert declinations and inclinations of two vectors (denoted as Dec1, Inc1, Dec2, Inc2) and click on the CALCULATE button (or press ENTER).
4. The orientation of the plane is displayed as Declination (Dec) and Inclination (Inc) in red text boxes; the Angle between vectors is in yellow text box. The plane together with the original lineations is displayed in the spherical projection as red circles and green and blue points, respectively.
5. To clear input boxes click on the CLEAR button (or press ESC).

5.3. Group statistic window

Group statistics window was designed for analyses on a group of specimens. It consists of:



- (1) Headline
- (2) Demagnetization States selector
- (3) Spherical projection plot
- (4) Coordinates selector
- (5) Display selector
- (6) Statistics selector
- (7) Records list
- (8) DISREGARD RECORD(S) button
- (9) RESTORE button
- (10) NORMAL button
- (11) REVERSED button
- (12) INVERT button
- (13) Mean direction frame
- (14) Orientation matrix frame
- (15) Virtual geomagnetic pole position (VGP) frame

◆ **Launching the Group statistics window**


Press *File* → *Group statistics* (or *CTRL + G*) in the main menu.

◆ **Closing the Group statistics window**

Press *File* → *Exit* (or *CTRL + X*) in the main menu.

◆ **Opening specimen files for the groups statistics**

1. Press *File* → *Open* (or *CTRL + O*) in the main menu. Open multiple window appears on screen.



2. Select type of file in the List file of types. Individual Specimen file(s), or Fisher statistics output file(s) (* .rsf) can be opened.
3. Select file(s) by either dragging the mouse over File names (left mouse button pressed [+CTRL for individual multiple selection, +SHIFT for continuous multiple selection]) or by clicking on the SELECT ALL button. To deselect all files click on the Clear selection button.
4. When desired files are selected click on the OPEN button.
5. After opening, the selected files are pooled together and all available demagnetization steps and principal component/great circle names are displayed in the States checkbox list.

◆ **Closing specimen files**

Press *File* → *Close* (or *CTRL + Z*) in the main menu. All selected files are closed and the States checkbox list is cleared.

◆ **Selecting demagnetization step/components/great circles**

1. Click on the States checkbox list to select demagnetization step(s) or component/great circle name(s) to be included in the statistics. More than one item can be selected (i.e. several differently-named components can be processed together).
2. Vector directions (denoted as Data in the Display frame) are displayed as black symbols, great circles as green arches. Mean direction (Mean) calculated from the currently displayed data and projection of the 95% confidence cone (A95) is plotted as magenta crossed symbol and magenta line, respectively. Solid symbols and lines represent directions in the lower hemisphere, open symbols and dashed lines represent directions in the upper hemisphere. Optionally, the projection of the semi-angle of precision/MAD of individual data (a95) can be displayed in gray line, and three eigenvectors of the orientation matrix (Orient. matrix) calculated from the currently displayed data plotted as blue numbers 1, 2, 3. Coordinate system corresponds to the highest system available in the selected files. (Note: Data for all coordinate systems may not be available in all selected files; the number of records may, thus, differ when different coordinate systems are chosen using the Coordinates selector).
3. List of data containing the specimen names, demagnetizing state, declination and inclination in the current coordinate system and precision/MAD is show in the Records list. Numerical results of the statistics calculated from currently displayed data are shown in the frames below the spherical projection.

Mean direction frame

- N – number of currently displayed data included in the statistics.
 - Dec – declination of the mean vector.
 - Inc – inclination of the mean vector.
 - R – length of the resultant vector (after vector addition of a set of *N* unit vectors); *R* is always $\leq N$ and approaches *N* only when the vectors are tightly clustered.
 - k – the best estimate of the precision parameter.
- $$k = \frac{N-1}{N-R}$$
- a95 – the angle within which the unknown true mean lies at 95% confidence level.

$$\cos \alpha_{(1-p)} = 1 - \frac{N-R}{R} \left\{ \left(\frac{1}{p} \right)^{\frac{1}{N-1}} - 1 \right\}$$

where $(1-p) = 0.95$, i.e. 95%.

◆ **Saving mean direction results**

1. Insert Name of the site and State (i.e. component name).
2. Click on the **SAVE** button.
3. Select/Insert a name of the Fisher statistics output file (*.rsf); when the file already exists the new record is appended in the end of the file.

Orientation matrix frame

- E1 – 1st eigenvalue and declination / inclination of the 1st eigenvector of orientation matrix.
- E2 – 2nd eigenvalue and declination / inclination of the 2nd eigenvector of orientation matrix.
- E3 – 3rd eigenvalue and declination / inclination of the 3rd eigenvector of orientation matrix.

Orientation matrix is calculated as:

$$O = \frac{1}{N} \begin{vmatrix} \sum l_i^2 & \sum l_i m_i & \sum l_i n_i \\ \sum l_i m_i & \sum m_i^2 & \sum m_i n_i \\ \sum l_i n_i & \sum m_i n_i & \sum n_i^2 \end{vmatrix}$$

where l_i , m_i , and n_i are the directional cosines of the i^{th} linear element.

◆ Saving the 3rd eigenvector results

1. Insert Name of the site and State (i.e. component/great circle name).
2. Click on the SAVE button.
3. Select/Insert a name of the Fisher statistics output file (*.rsf); when the file already exist the new record is appended in the end of the file.

♥**Tip:** The 3rd eigenvector of the orientation matrix can be used as an estimate of the intersection of a number of great circles.

VGP (Virtual geomagnetic pole) frame

- Site Lat – geographic latitude of the sampling site. Default location is set in the Setting window; after inserting a new location press the Recalculate button. New location is automatically set as default site location.
- Site Lon – geographic longitude of the sampling site
Sign conventions for geographic locations: Latitudes increase from -90° at south geographic pole to 0° at equator and to $+90^\circ$ at the north geographic pole. Longitudes east of the Greenwich meridian are positive, while westerly longitudes are negative.
- Pole Lat – geographical latitude of the virtual pole
- Pole Lon – geographical longitude of the virtual pole
- PaleoLat – magnetic colatitude which is the great circle distance from site to pole
- dp – semi-axis of the confidence ellipse along the great circle path from site to pole
- dm – semi-axis of the confidence ellipse perpendicular to that great circle path
- ASD – angular standard deviation of individual VGPs from the mean VGP

◆ Saving VGP results

1. Insert Name of the site
2. Click on the SAVE button
3. Select/Insert a name of the VGPs output file (*.txt); when the file already exists the new record is appended in the end of the file.

◆ **Subsequent selecting within displayed data**

1. Drag a rectangle in the spherical projection plot: left mouse button pressed selects all data points within the rectangle, right mouse button pressed selects data points projected on the lower hemisphere, middle mouse button pressed selects data points projected on the upper hemisphere. Selected data points appear in red color in the plot and are also selected in the Records list. To deselect all data points click anywhere in the plot.
2. Click the mouse on the desired record in the Records list. For multiple selection drag mouse with the left mouse button pressed over desired records in the Records list [+CTRL for individual multiple selection, +SHIFT for continuous multiple selection]. Selected records appear as the red data point in the spherical projection plot. To deselect all records click on the headline of the Records list.
3. To invert selection click on the INVERT button.
4. After calculating the orientation matrix the whole displayed data set can be divided into two groups of data. One groups contains vectors deflected less than 90° from the first eigenvector of the orientation matrix while second groups contains vectors deflected more than 90° from the first eigenvector of the orientation matrix. The groups are technically denoted as normally- polarized data and reversely-polarized data, respectively. To select the normally-polarized data click on the NORMAL button. To select the reversely-polarized data click on the REVERSED button.

◆ **Disregarding selected record**

1. Selected record(s) can be disregarded from the statistics by clicking on the DISREGARD RECORD(S) button. All group statistics results are automatically recalculated (i.e. Mean direction, Orientation matrix, and VGP). Multiple disregarding can be done. Note that the record are disregards in the virtual memory file only; no records are deleted from the specimen files store on disk!
2. Click on the RESTORE button to restore the original set of records.

◆ **Transposition of the reversely-polarized data**

By default both normally- and reversely-polarized vectors are displayed in their actual orientation. Click on the Transposed option in the Statistic frame in order to display the reversed vectors of the reversely-polarized data. Click on the Normal option to display all data in their actual orientation. Note: After transposition all group statistics results are automatically recalculated (i.e. Mean direction, and VGP).

♥ **Tip:** Group statistics is calculated from the currently displayed records. Combined use of the mouse selections; Normal, Reversed, Invert, Disregard record(s), and Restore buttons; and Normal and Transposed statistics options provides a power tool for calculating group statistics from various data subsets. Statistics results of different data subsets can be appended into the output files using the Save buttons; currently displayed graphics can be either printed or exported to the graphic file.

◆ **Saving individual data**

1. Press *File* → *Save as* (or *CTRL + S*) in the main menu.
2. Select type of file: Group statistics file (*.rsf), Magnetostratigraphy MPS file (*.txt), SpheriStat file (*.dat), or original Agico JR6 file (*.jr6).


3. Insert file name and click on the SAVE button.

◆ **Exporting the graphics into graphic file**

1. Press *Graphics* → *Export* (or *CTRL + E*) in the main menu.
2. Select directory and file name.
3. Graphics will be stored as a vector graphics in the Windows metafile allowing for further image processing.

◆ **Printing the currently displayed graphics**

1. Press *Graphics* → *Print* (or *CTRL + P*) in the main menu.
2. Select a printer
Note: Selected printer will be set as Windows default printer.
3. Printing graphics using: window appears. Specify printing options.




4. Click on the PRINT button. To abort printing click on the CANCEL button.

5.3.1. Tool windows

5.3.1.1. Angle between vectors

◆ **Calculating angle between two vectors**

1. Press *Tools* → *Angle between vectors* in the main menu.
2. Angle between vectors window appears.




3. Manually insert declinations and inclinations of two vectors (denoted as Dec1, Inc1, Dec2, Inc2) and click on the CALCULATE button.
4. Result is displayed in the green box (denoted as Angle).
5. To clear input boxes click on the CLEAR button.

5.3.1.2. Magnetic inclination

◆ Calculating magnetic inclination corresponding to geographic latitude and vice versa

1. Press *Tools* → *Magnetic inclination* in the main menu.
2. Magnetic inclination window appears.




3. Manually insert either geographic Latitude or magnetic Inclination and click on the appropriate arrow.
4. To clear input boxes click on the CLEAR button.

REFERENCES


- Fisher, R. A., 1953. Dispersion on a sphere. *Proc. R. Soc. London A*, **127**, 295–305.
- Halls, H. C., 1976. A least-squares method to find a remanence direction from converging remagnetization circles. *Geophys. J. R. Astr. Soc.*, **45**, 297–304.
- Kirschvink, J. L., 1980. The least-squares line and plane and the analyses of paleomagnetic data. *Geophys. J. R. Astr. Soc.*, **62**, 699–718.
- Lowrie, W., 1990. Identification of ferromagnetic minerals in a rock by coercivity and unblocking temperature properties. *Geoph. Res. Lett.*, **2**, 159–162.
- Zijderveld, J. D. A., 1967. A.C. demagnetisation methods. In: Runcorn, S.D., Collinson, D. W. (eds.), *Methods in Palaeomagnetism*. Elsevier, Amsterdam, 254–286.

ORIENTATION PARAMETERS

Parameter P1



It is clock value of **the direction** of the fiducial **mark drawn** on the frontal side of cylinder. This arrow is X1 axis of the specimen coordinate system. The orientation of the arrow may, or need not, be measured.



P1=12

P1=3

P1=6

P1=9

Parameter P2

Its value is **0** or **90**.

P2=0 if the dip of the frontal side is measured.

P2=90 if the plunge of the cylinder (drilling) axis is measured.

Parameter P3

It is clock value of **the direction** (visualized by arrow which need not necessarily be drawn) **which is measured** in the field.

Parameter P4

P4=0 Value zero means that azimuth of dip and dip of mesoscopic foliation are measured.

P4=90 Value 90 means that strike (right oriented) and dip are measured.

Examples

Agico system:

P1 = 12	fiducial mark is oriented upwards
P2 = 90	plunge of cylinder axis is measured
P3 = 6	azimuth of dip of frontal plane is measured in the field
P4 = 0	azimuth of dip and dip of mesoscopic foliation are measured

University of Santa Barbara:

P1 = 12	fiducial mark is oriented upwards
P2 = 0	dip of frontal line is measured
P3 = 3	the strike of frontal plane is measured in the field
P4 = 90	strike and dip of mesoscopic foliation are measured

Paleomagnetic laboratory in Espoo:

P1 = 12	fiducial mark is oriented upwards and its orientation is also measured in the field => P1 = P3 = 12
P2 = 90	plunge of cylinder axis is measured
P3 = 12	see P1
P4 = 0	azimuth of dip and dip of mesoscopic foliation are measured